


# K3k

# THREE FUNCTIONS MANUAL BED

Registration code: SK-CD1



**Healthward** uk

Healthier your life

# K3k

# FEATURES AND HIGHLIGHTS









# Technical parameters

External size (LxW):	2105x980 mm
Mattress platform:	1925x820 mm
Safe working load:	200 kg
Height adjustment:	400-710mm
Back-rest adjustment:	0-75° (±10°)
Knee-rest adjustment:	0-40° (±10°)

# Technical configuration

☑ Manual crank	3 pcs
☑ 5" covered castors	4 pcs
☑ Foldable side rail	1 set
☑ ABS Bed Ends with Safe Lock	1set
☑ IV Pole	1set
✓ IV Pole Prevision	4pcs
☑ Drainage Hook	4pcs
✓ Mattress Holder	2pcs

# Optional color (MOQ100pcs-color changing)



Red oak





#### Design

The design is both modern and practical.
The head and foot boards feature

a beautiful and elegant shape. The use of a hollow protective wrap angle, combined with a barreltype bumper, makes the product instantly recognizable.



#### Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



#### Plastic | Painting

Environmental plastic can be sterilized at a temperature of 100°C and has a tensile strength of up to

11 Process epoxy painting with a thickness of 0.12mm and a brightness of 60°, which is capable of resisting impacts of up to 50kg.

# SAIKANO

# **Accessories**



#### **Head & foot board**

The head and foot boards feature "integrated protective corners" with hollow design, which are combined with barrel-type bumper wheels that have a diameter of  $60\,\mathrm{mm}$ .

This design not only minimizes the overall space occupied by the bed but also provides protection.

Additionally, the head and foot boards are effortless to disassemble and clean



#### Castors

5" covered casters.

PR material, walks 30KM with no wear and tear. Precision-designed antiwinding shell with no screws.

Pass Dynamic test: withstand a weight of 120kg, walk 30KM, and pass through 500 obstacles without any deformation or damage.















Back-rest adjustment

75°



Knee-rest adjustment

40°



Height adjustment

 $400 - 710_{\mathsf{mm}}$ 















Design

Research

Manufacturing

Maintenance

Service



#### USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

## United Kingdom

**UK Healthward International Ltd.**Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

#### Russia

**Saikang Medical Russia LLC** Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 2145566 / 8811 / 1516

### East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254733704690

### Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

# **Healthward** Tuk



- No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road · Wuhou District (Chenedu) · China
- +86 18021231901 · \$\sum\_2\$ 24/7
- export@health-ward.com
- www.health-ward.com







Company video

General catalogue